

iGenomX Sample Demultiplexing Tool

This document describes how to use a command line tool for sample demultiplexing iGenomX data and is intended for individuals familiar with using a command line bioinformatics tool.

Table of Contents

- [Example Command Line](#)
- [Prerequisites](#)
 - [Java 8](#)
 - [Fgbio](#)
- [Running the Tool](#)
 - [Inputs](#)
 - [Outputs](#)
 - [Advanced Options](#)
- [Appendix](#)
 - [Additional Recommended Best Practices](#)
 - [Converting BAMs to FASTQ](#)
 - [Memory Usage](#)
 - [Example Sample Sheet](#)
 - [Example Metadata CSV](#)
 - [Example Sample Barcode Metrics](#)

Example Command Line

Below shows an example command line for a paired end 2x150bp experiment:

```
java -Xmx8G -jar fgbio.jar DemuxFastqs \  
  --inputs example_S1_L001_R1_001.fastq.gz example_S1_L001_R2_001.fastq.gz \  
  --metadata SampleSheet.csv \  
  --read-structures 8B8M134T 8M142T \  
  --output /path/to/output/directory \  
  --metrics example.sample_barcode_metrics.txt
```

A BAM per sample will be written to the `path/to/output/directory` directory, including a BAM for reads that were not assigned a sample.

The `--output-fastqs true` option can be outputted to output FASTQs instead of BAMs.

Prerequisites

Java 8

[Java 8](#) is required to run the command line tools.

Fgbio

The `fgbio` tool set can be downloaded from <https://github.com/fulcrumgenomics/fgbio>. Version 0.2.0 or higher is required; at the time of this writing, building from source is required.

To install from source, see [Building Fgbio](#).

Please note that `git`, `java8`, the java JDK, `r-base`, `sbt`, and `scala` are required to build from source. The executable JAR can be found in `fgbio/target/scala-2.11/fgbio-*.jar` after building.

Fgbio is an open source tool developed by [Fulcrum Genomics](#). Please visit the [Fgbio project page](#) to obtain support and report any bugs.

Running the Tool

The `DemuxFastqs` tool will be used within the `fgbio` set of tools. Please see the [Prerequisites](#) section for obtaining the tool set. To view the available commands, use `java -jar fgbio.jar --help`. To view the all options for the `DemuxFastqs` tool, use `java -jar fgbio.jar DemuxFastqs --help`.

Inputs

Input FASTQs

The input FASTQ(s) should be provided with the `--inputs` option. For paired end reads, the FASTQ corresponding to read one and read two should be specified in sequence:

```
--inputs <read_one>.fastq.gz <read_two>.fastq.gz
```

Compressed FASTQs (gzipped: `fastq.gz`) are supported.

Read Structure

Metadata CSV

In lieu of a [Sample Sheet](#), a CSV file with just the sample section (`[Data]`) described in the [Sample Sheet](#) can be given.

An [Example Metadata CSV](#) is shown in the [Appendix](#).

Outputs

The output directory should be specified with the `--output` option. The output directory will contain one BAM file per sample in the sample sheet, plus a BAM for reads that could not be assigned to a sample given the criteria.

The BAM file name for a sample will be the concatenation of sample id, sample name, and sample barcode bases (from the sample sheet), delimited by “-” (i.e.

`<sample_id>-<sample_name>-<sample_barcode>.bam`). The BAM’s read group will have sample id, sample name, and library id corresponding to the similarly named values in the sample sheet. The library id will be the sample id if not found, and the platform unit will be the sample name concatenated with the sample barcode bases delimited by a “.” (i.e. `<sample_id>.<sample_barcode>.bam`). Additional command line options are available to specify additional metadata in the BAM’s read group. The name for the unmatched sample is `unmatched` by default, but can be specified using the `--unmatched` option.

Alternatively, gzipped FASTQs can be written using the “`-output-fastqs=true`” option instead of BAMs.

For paired end data, the output will have the suffix “`R1.fastq.gz`” and “`R2.fastq.gz`” for read one and read two respectively. The sample barcode and molecular barcodes (concatenated) will be appended to the read name and delimited by a colon.

A metrics file will also be output providing analogous information to the metric described here:

[SampleBarcodeMetric](#). Use the `--metrics` option to specify the path to where the metrics file should be written. An `[Example Sample Barcode Metrics][#example-sample-barcode-metrics]` file is show in the [Appendix](#).

Advanced Options

For a read to be assigned to a sample two criteria must be met:

1. The read’s sample barcode must match the sample’s barcode with `<= the maximum allowed number of mismatches` .
2. The read must not match any other sample’s barcode with `<= mismatches to best barcode + minimum mismatch delta` .

For example, with `--max-mismatches=1` and `--min-delta=2` :

- If a sample is matched with a perfect match (no mismatches) all other samples' barcodes must be at least 2 mismatches away from the barcode sequence.
- If a sample is matched with a single mismatch, all other samples' barcodes must be at least 3 mismatches away from the barcode sequence.

Additional options are shown in the usage by specifying the `--help` option.

Appendix

Additional Recommended Best Practices

Adapter Marking

Care must be taken when adapter trimming or adapter marking. If the insert size (length of fragment being sequenced) is less than the number of bases sequence, both the inline sample barcode bases for the opposite read pair and the adapter sequenced may be present at the end of the given read. Alignment tools like `bwa` will soft-clip bases at the ends of reads that do not map well. Post-alignment tools like `MergeBamAlignment` (`Picard`) or `ClipOverlappingReads` (`fgbio`).

Specifying the adapter sequences to match, for example with `Picard` tool `MarkIlluminaAdapters` to include leading masked bases (`N` s) may be appropriate. Furthermore, using an alignment tool that can soft-clip the ends of reads (ex. `bwa`) can also be used.

PCR Duplicate Marking

It is recommended to use the molecular barcodes for each read when marking PCR duplicates, for example with the `Picard` tool `MarkDuplicates` by specifying the `BARCODE_TAG=RX` command line option. This allow duplicate marking to better discriminate true PCR duplicates versus those that mapped to the same location by chance, which can occur frequently when sequencing to deep coverage.

Converting BAMs to FASTQ

In some cases it is necessary to convert a BAM to FASTQ for downstream processing. The `Picard` tool `SamToFastq` is recommended for this conversion. Since the sample barcode AND molecular barcodes are stored in SAM tags, they will not be present in the converted FASTQs using this tool. This, if FASTQs are required for mapping (ex. `bwa`), it is recommended to use the `Picard` tool `MergeBamAlignment` to post-

process the mapped BAM and thus restore the various read-level metadata (ex. SAM tags and read groups) stored in the original demultiplexed BAM.

Alternatively, the “`--output-fastqs=true`” option can be used to write gzipped FASTQs instead of BAMs, but with the loss of sample metadata that can be stored in a read group in a BAM’s header.

Memory Usage

It is recommended to run the tool with at least 8GB of memory, which can be specified when running the tool:

```
java -Xmx8g -jar fgbio.jar ... .
```

Example Sample Sheet

Below is an minimal example sample sheet for demultiplexing 96 samples.

```
[Header],,,,,,,
IEMFileVersion,4,,,,,,
Investigator Name,Joe,,,,,,
Experiment Name,EXPID,,,,,,
Date,1/1/00,,,,,,
Workflow,GenerateFASTQ,,,,,,
Application,FASTQ Only,,,,,,
Assay,Assay Name,,,,,,
Description,The Description,,,,,,
Chemistry,Amplicon,,,,,,

,,,,,,
[Reads],,,,,,,
151,,,,,,
151,,,,,,

,,,,,,
[Settings],,,,,,,
ReverseComplement,0,,,,,,
Adapter,AGATCGGAAGAGCACACGTCTGAACTCCAGTCA,,,,,,
AdapterRead2,AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT,,,,,,

,,,,,,
[Data],,,,,,,
Sample_ID,Sample_Name,Library_Id,Sample_Plate,Sample_Well,index,Sample_Project,Descri
ption
ID-1,1,,,,,AACCAAGG,,
ID-2,2,,,,,AACGTTGC,,
ID-3,3,,,,,AAGGTAGC,,
ID-4,4,,,,,ACACGTGT,,
ID-5,5,,,,,ACCAAGGA,,
```

ID-6,6,,,,ACCTGTTC,,
ID-7,7,,,,ACGTAGGA,,
ID-8,8,,,,ACTCGTGA,,
ID-9,9,,,,AGACAGTG,,
ID-10,10,,,,AGAGGTGT,,
ID-11,11,,,,AGCTAGGA,,
ID-12,12,,,,AGTCAGTC,,
ID-13,13,,,,AGTGGTCT,,
ID-14,14,,,,AGTGGTGA,,
ID-15,15,,,,ATCCTACC,,
ID-16,16,,,,ATCCTAGG,,
ID-17,17,,,,ATGCGCAT,,
ID-18,18,,,,ATGCGCTA,,
ID-19,19,,,,ATTAGCGC,,
ID-20,20,,,,ATTAGGCC,,
ID-21,21,,,,CAACGTTG,,
ID-22,22,,,,CAAGCTAC,,
ID-23,23,,,,CAAGCTTG,,
ID-24,24,,,,CACATGAC,,
ID-25,25,,,,CACATGTG,,
ID-26,26,,,,CAGACTCA,,
ID-27,27,,,,CAGACTGT,,
ID-28,28,,,,CAGTTGAC,,
ID-29,29,,,,CAGTTGTG,,
ID-30,30,,,,CATGCTAG,,
ID-31,31,,,,CATGCTTC,,
ID-32,32,,,,CCATATCC,,
ID-33,33,,,,CCATATGG,,
ID-34,34,,,,CCTACCAT,,
ID-35,35,,,,CCTACCTA,,
ID-36,36,,,,CGAACCAT,,
ID-37,37,,,,CGAACCTA,,
ID-38,38,,,,CGATTACG,,
ID-39,39,,,,CGATTAGC,,
ID-40,40,,,,CGTAGCAT,,
ID-41,41,,,,CGTAGCTA,,
ID-42,42,,,,CTACCAAG,,
ID-43,43,,,,CTACCATC,,
ID-44,44,,,,CTAGTCCT,,
ID-45,45,,,,CTAGTCGA,,
ID-46,46,,,,CTCTCACA,,
ID-47,47,,,,CTCTCAGT,,
ID-48,48,,,,CTGATCAG,,
ID-49,49,,,,CTGATCTC,,
ID-50,50,,,,CTTCCAAC,,

ID-51,51,,,,CTTCCATG,,
ID-52,52,,,,CTTGTCCA,,
ID-53,53,,,,CTTGT CGT,,
ID-54,54,,,,GAAGCAAC,,
ID-55,55,,,,GAAGCATG,,
ID-56,56,,,,GACATCAC,,
ID-57,57,,,,GACATCTG,,
ID-58,58,,,,GAGACACA,,
ID-59,59,,,,GAGACAGT,,
ID-60,60,,,,GAGTTCAC,,
ID-61,61,,,,GAGTTCTG,,
ID-62,62,,,,GATGCAAG,,
ID-63,63,,,,GATGCATC,,
ID-64,64,,,,GCATAACC,,
ID-65,65,,,,GCATAAGG,,
ID-66,66,,,,GCGCTATA,,
ID-67,67,,,,GCGCTTAA,,
ID-68,68,,,,GCTTCCAT,,
ID-69,69,,,,GCTTCCTA,,
ID-70,70,,,,GGATCCAT,,
ID-71,71,,,,GGATCCTA,,
ID-72,72,,,,GGTACGAA,,
ID-73,73,,,,GGTACGTT,,
ID-74,74,,,,GTACAGCT,,
ID-75,75,,,,GTACAGGA,,
ID-76,76,,,,GTAGGTAG,,
ID-77,77,,,,GTCTAGAC,,
ID-78,78,,,,GTGAGTCT,,
ID-79,79,,,,G TTCAGCA,,
ID-80,80,,,,GTTGGTAC,,
ID-81,81,,,,TACGAACC,,
ID-82,82,,,,TAGCTACG,,
ID-83,83,,,,TATTGCGG,,
ID-84,84,,,,TCAGCTCT,,
ID-85,85,,,,TCCATGCT,,
ID-86,86,,,,TCGACTAG,,
ID-87,87,,,,TCGTTGCT,,
ID-88,88,,,,TCTGCTCA,,
ID-89,89,,,,T GACTGAC,,
ID-90,90,,,,TGCACTAG,,
ID-91,91,,,,TGCTTGCT,,
ID-92,92,,,,TGGTCTAG,,
ID-93,93,,,,TGTCTGAG,,
ID-94,94,,,,TTATCCGC,,
ID-95,95,,,,TTCGGCAA,,

```
ID-96,96,,,,TTGGCCAA,,
```

Example Metadata CSV

Below is an minimal example metadta CSV file for demultiplexing 96 samples.

```
Sample_ID,Sample_Name,Library_Id,Sample_Plate,Sample_Well,index,Sample_Project,Description
ID-1,1,,,,AACCAAGG,,
ID-2,2,,,,AACGTTGC,,
ID-3,3,,,,AAGGTAGC,,
ID-4,4,,,,ACACGTGT,,
ID-5,5,,,,ACCAAGGA,,
ID-6,6,,,,ACCTGTTC,,
ID-7,7,,,,ACGTAGGA,,
ID-8,8,,,,ACTCGTGA,,
ID-9,9,,,,AGACAGTG,,
ID-10,10,,,,AGAGGTGT,,
ID-11,11,,,,AGCTAGGA,,
ID-12,12,,,,AGTCAGTC,,
ID-13,13,,,,AGTGGTCT,,
ID-14,14,,,,AGTGGTGA,,
ID-15,15,,,,ATCCTACC,,
ID-16,16,,,,ATCCTAGG,,
ID-17,17,,,,ATGCGCAT,,
ID-18,18,,,,ATGCGCTA,,
ID-19,19,,,,ATTAGCGC,,
ID-20,20,,,,ATTAGGCC,,
ID-21,21,,,,CAACGTTG,,
ID-22,22,,,,CAAGCTAC,,
ID-23,23,,,,CAAGCTTG,,
ID-24,24,,,,CACATGAC,,
ID-25,25,,,,CACATGTG,,
ID-26,26,,,,CAGACTCA,,
ID-27,27,,,,CAGACTGT,,
ID-28,28,,,,CAGTTGAC,,
ID-29,29,,,,CAGTTGTG,,
ID-30,30,,,,CATGCTAG,,
ID-31,31,,,,CATGCTTC,,
ID-32,32,,,,CCATATCC,,
ID-33,33,,,,CCATATGG,,
ID-34,34,,,,CCTACCAT,,
ID-35,35,,,,CCTACCTA,,
ID-36,36,,,,CGAACCAT,,
```

ID-37,37,,,,CGAACCTA,,
ID-38,38,,,,CGATTACG,,
ID-39,39,,,,CGATTAGC,,
ID-40,40,,,,CGTAGCAT,,
ID-41,41,,,,CGTAGCTA,,
ID-42,42,,,,CTACCAAG,,
ID-43,43,,,,CTACCATC,,
ID-44,44,,,,CTAGTCCT,,
ID-45,45,,,,CTAGTCGA,,
ID-46,46,,,,CTCTCACA,,
ID-47,47,,,,CTCTCAGT,,
ID-48,48,,,,CTGATCAG,,
ID-49,49,,,,CTGATCTC,,
ID-50,50,,,,CTTCCAAC,,
ID-51,51,,,,CTTCCATG,,
ID-52,52,,,,CTTGTCCA,,
ID-53,53,,,,CTTGTCGT,,
ID-54,54,,,,GAAGCAAC,,
ID-55,55,,,,GAAGCATG,,
ID-56,56,,,,GACATCAC,,
ID-57,57,,,,GACATCTG,,
ID-58,58,,,,GAGACACA,,
ID-59,59,,,,GAGACAGT,,
ID-60,60,,,,GAGTTCAC,,
ID-61,61,,,,GAGTTCTG,,
ID-62,62,,,,GATGCAAG,,
ID-63,63,,,,GATGCATC,,
ID-64,64,,,,GCATAACC,,
ID-65,65,,,,GCATAAGG,,
ID-66,66,,,,GCGCTATA,,
ID-67,67,,,,GCGCTTAA,,
ID-68,68,,,,GCTTCCAT,,
ID-69,69,,,,GCTTCCTA,,
ID-70,70,,,,GGATCCAT,,
ID-71,71,,,,GGATCCTA,,
ID-72,72,,,,GGTACGAA,,
ID-73,73,,,,GGTACGTT,,
ID-74,74,,,,GTACAGCT,,
ID-75,75,,,,GTACAGGA,,
ID-76,76,,,,GTAGGTAG,,
ID-77,77,,,,GTCTAGAC,,
ID-78,78,,,,GTGAGTCT,,
ID-79,79,,,,GTTTCAGCA,,
ID-80,80,,,,GTTGGTAC,,
ID-81,81,,,,TACGAACC,,

6	6	ACCTGTTC	11772	11772	11586	11586	
	186	186			0.000479	0.011847	
46418		0.000479	0.011847				0.0
7	7	ACGTAGGA	9182	9182	9107	9107	
	75	75			0.000373	0.00924	
36205		0.000373	0.00924				0.0
8	8	ACTCGTGA	9446	9446	9390	9390	
	56	56			0.000384	0.009506	
37246		0.000384	0.009506				0.0
9	9	AGACAGTG	709179	709179	702664	702664	
	6515	6515			0.028828	0.713675	
96338		0.028828	0.713675				2.7
10	10	AGAGGTGT	596421	596421	593342	593342	
	3079	3079			0.024245	0.600202	
51726		0.024245	0.600202				2.3
11	11	AGCTAGGA	351888	351888	346499	346499	
	5389	5389			0.014304	0.354119	
87517		0.014304	0.354119				1.3
12	12	AGTCAGTC	93779	93779	92942	92942	
	837	837			0.003812	0.094374	
69777		0.003812	0.094374				0.3
13	13	AGTGGTCT	368605	368605	365969	365969	
	2636	2636			0.014984	0.370942	
53433		0.014984	0.370942				1.4
14	14	AGTGGTGA	187000	187000	185816	185816	
	1184	1184			0.007602	0.188186	
37353		0.007602	0.188186				0.7
15	15	ATCCTACC	592621	592621	588125	588125	
	4496	4496			0.02409	0.596378	
36742		0.02409	0.596378				2.3
16	16	ATCCTAGG	158134	158134	156775	156775	
	1359	1359			0.006428	0.159137	
23532		0.006428	0.159137				0.6
17	17	ATGCGCAT	405648	405648	401566	401566	

	4082		4082			0.01649	0.40822	
						0.01649	0.40822	1.5
99496								
18	18		ATGCGCTA	329727	329727	326387	326387	
	3340		3340			0.013403	0.331817	
						0.013403	0.331817	1.3
00134								
19	19		ATTAGCGC	993700	993700	986891	986891	
	6809		6809			0.040394	1	
						0.040394	1	3.9
18222								
20	20		ATTAGGCC	21223	21223	20057	20057	
	1166		1166			0.000863	0.021358	
						0.000863	0.021358	0.0
83684								
21	21		CAACGTTG	379670	379670	376922	376922	
	2748		2748			0.015434	0.382077	
						0.015434	0.382077	1.4
97063								
22	22		CAAGCTAC	373005	373005	371689	371689	
	1316		1316			0.015163	0.37537	
						0.015163	0.37537	1.4
70782								
23	23		CAAGCTTG	10444	10444	10425	10425	
	19		19			0.000425	0.01051	
						0.000425	0.01051	0.0
41181								
24	24		CACATGAC	338900	338900	338429	338429	
	471		471			0.013776	0.341049	
						0.013776	0.341049	1.3
36304								
25	25		CACATGTG	13034	13034	13003	13003	
	31		31			0.00053	0.013117	
						0.00053	0.013117	0.0
51394								
26	26		CAGACTCA	250271	250271	249225	249225	
	1046		1046			0.010174	0.251858	
						0.010174	0.251858	0.9
86834								
27	27		CAGACTGT	7349	7349	7314	7314	
	35		35			0.000299	0.007396	
						0.000299	0.007396	0.0
28978								
28	28		CAGTTGAC	155718	155718	155431	155431	
	287		287			0.00633	0.156705	

51	51	CTTCCATG	9826	9826	9566	9566	
	260	260			0.000399	0.009888	
			0.000399	0.009888			0.0
38745							
52	52	CTTGTCCTA	235932	235932	233301	233301	
	2631	2631			0.009591	0.237428	
			0.009591	0.237428			0.9
30295							
53	53	CTTGTCGT	212940	212940	211585	211585	
	1355	1355			0.008656	0.21429	
			0.008656	0.21429			0.8
39636							
54	54	GAAGCAAC	17216	17216	17146	17146	
	70	70			0.0007	0.017325	
			0.0007	0.017325			0.0
67884							
55	55	GAAGCATG	725654	725654	722869	722869	
	2785	2785			0.029498	0.730255	
			0.029498	0.730255			2.8
613							
56	56	GACATCAC	17280	17280	17192	17192	
	88	88			0.000702	0.01739	
			0.000702	0.01739			0.0
68136							
57	57	GACATCTG	599177	599177	598211	598211	
	966	966			0.024357	0.602976	
			0.024357	0.602976			2.3
62593							
58	58	GAGACACA	16878	16878	16783	16783	
	95	95			0.000686	0.016985	
			0.000686	0.016985			0.0
66551							
59	59	GAGACAGT	588025	588025	585325	585325	
	2700	2700			0.023903	0.591753	
			0.023903	0.591753			2.3
1862							
60	60	GAGTTCAC	16337	16337	16237	16237	
	100	100			0.000664	0.016441	
			0.000664	0.016441			0.0
64418							
61	61	GAGTTCTG	176093	176093	175633	175633	
	460	460			0.007158	0.177209	
			0.007158	0.177209			0.6
94346							
62	62	GATGCAAG	741904	741904	738970	738970	

	2934		2934			0.030158	0.746608	
						0.030158	0.746608	2.9
25374								
63	63		GATGCATC	138194	138194	137521	137521	
	673		673			0.005618	0.13907	
						0.005618	0.13907	0.5
44908								
64	64		GCATAACC	708525	708525	706272	706272	
	2253		2253			0.028802	0.713017	
						0.028802	0.713017	2.7
93759								
65	65		GCATAAGG	354180	354180	352832	352832	
	1348		1348			0.014397	0.356425	
						0.014397	0.356425	1.3
96554								
66	66		GCGCTATA	475894	475894	470844	470844	
	5050		5050			0.019345	0.478911	
						0.019345	0.478911	1.8
7648								
67	67		GCGCTTAA	453613	453613	448494	448494	
	5119		5119			0.018439	0.456489	
						0.018439	0.456489	1.7
88625								
68	68		GCTTCCAT	377381	377381	376226	376226	
	1155		1155			0.015341	0.379774	
						0.015341	0.379774	1.4
88037								
69	69		GCTTCCTA	9025	9025	8851	8851	
	174		174			0.000367	0.009082	
						0.000367	0.009082	0.0
35586								
70	70		GGATCCAT	129093	129093	128657	128657	
	436		436			0.005248	0.129911	
						0.005248	0.129911	0.5
09022								
71	71		GGATCCTA	7676	7676	7573	7573	
	103		103			0.000312	0.007725	
						0.000312	0.007725	0.0
30267								
72	72		GGTACGAA	124991	124991	123288	123288	
	1703		1703			0.005081	0.125783	
						0.005081	0.125783	0.4
92847								
73	73		GGTACGTT	8126	8126	8029	8029	
	97		97			0.00033	0.008178	

			0.00033		0.008178			0.0
32041								
74	74		GTACAGCT	211000	211000	209308	209308	
	1692		1692			0.008577	0.212338	
						0.008577	0.212338	0.8
31986								
75	75		GTACAGGA	14573	14573	14419	14419	
	154		154			0.000592	0.014665	
						0.000592	0.014665	0.0
57462								
76	76		GTAGGTAG	99705	99705	98679	98679	
	1026		1026			0.004053	0.100337	
						0.004053	0.100337	0.3
93143								
77	77		GTCTAGAC	10939	10939	10816	10816	
	123		123			0.000445	0.011008	
						0.000445	0.011008	0.0
43133								
78	78		GTGAGTCT	9840	9840	9678	9678	
	162		162			0.0004	0.009902	
						0.0004	0.009902	0.0
388								
79	79		G TTCAGCA	16317	16317	16164	16164	
	153		153			0.000663	0.01642	
						0.000663	0.01642	0.0
64339								
80	80		G TTGGTAC	9392	9392	9278	9278	
	114		114			0.000382	0.009452	
						0.000382	0.009452	0.0
37033								
81	81		TACGAACC	640623	640623	634304	634304	
	6319		6319			0.026041	0.644685	
						0.026041	0.644685	2.5
26017								
82	82		TAGCTACG	509337	509337	504551	504551	
	4786		4786			0.020705	0.512566	
						0.020705	0.512566	2.0
08348								
83	83		TATTGCGG	636050	636050	629482	629482	
	6568		6568			0.025856	0.640083	
						0.025856	0.640083	2.5
07985								
84	84		TCAGCTCT	286300	286300	283936	283936	
	2364		2364			0.011638	0.288115	
						0.011638	0.288115	1.1

96	96	TTGGCCAA	359032	359032	354958	354958
	4074	4074			0.014595	0.361308
		0.014595	0.361308			1.4
15686						
unmatched	unmatched	NNNNNNNN	525880	525880	0	0
	0	0			0.021377	0.529214
		0.021377	0.529214			2.0
73578						